

# ONESPRAY ALCOHOL BASED FAST ACTING SPRAY DISINFECTANT FOR MEDICAL EQUIPMENTS



EN 13624  
EN 13727  
EN 14476  
EN 14348

## USAGE

OneSpray is an alcohol-based product that is used for rapid disinfection and cleaning on medical device surfaces. It is ready to use. It is sprayed to cover the area to be disinfected and its entire surface from 30cm away. Should be waited for antimicrobial activity to occur, then wiped with a clean cloth/cotton. It does not contain aldehyde and phenol. It is based on ethanol, propanol and quaternary ammonium combination and does not fade or dye the material used. It has a nice fragrance.

## PROPERTIES

Ready to use product  
Aldehyde and phenol free  
Does not do color fading or dyeing on the material used  
Pleasant-smelling



**%99.9**  
Effective



## MICROBIOLOGICAL EFFECT

Bactericidal, Fungicidal, and Virucidal effects.

## COMPOSITION

%30 Ethanol, %10 2-Propanol, %0,25 Didecylmethylpoly(oxethyl) Ammonium Propionate, Auxiliary Substances, Essences, Deionized Water

## WARNINGS AND STORAGE CONDITIONS

Keep out of reach of children.  
Keep away from sunlight and high temperatures.  
Store in tightly closed original packaging/container, dry and cool area.  
Keep away from heat, sparks and fire.  
Under normal conditions, shelf life is 2 years.  
Use opened product within 6 months.  
Store locked up.



Alcohol Based



Ready to use



Aldehyde and phenol free



Pleasant-smelling



Fast and effective disinfection

KAF G20 - 1	750 ML
KAF G20 - 2	1000 ML
KAF G20 - 3	5 L

	Test Method	Microorganism	Effect Duration
<b>Bactericidal</b>	EN 13727	Staphylococcus aureus	<b>30 sec</b>
		Pseudomonas aeruginosa	
		Enterococcus hirae	
<b>Fungicidal / Yeas-ticidal</b>	EN 13624	Candida albicans	<b>5 min</b>
		Aspergillus niger	
<b>Tuberculocidal / Mycobactericidal</b>	EN 14348	Mycobacterium avium	<b>1 min</b>
		Mycobacterium terrae	
<b>Virucidal</b>	EN 14476	Poliovirus Type 1	<b>1 min</b>
		Human Adenovirus Type 5	
		Murine Norovirus	