

# OneGel

## LUBRICANT GEL WITH LIDOCAINE (STERILE)

### USAGE

Sterile Lubricating gel with lidocaine for mucosal instillation.  
Sterile gel is used for lubricating urethral catheters and other medical devices.  
during urethral application of therapeutic devices such as permanent catheters or  
dilatation, diagnostic devices such as ureteroscopy, endoscopy or cystoscopy.

### PROPERTIES

Sterile and disposable product.  
Has antiseptic properties.  
Has anesthetic effect.  
Eliminates the risk of damage to tissue.

### USAGE

Carefully insert the cone of the syringe, and slowly instill the sterile lubricant  
gel.



12,5 gr



### COMPOSITION

Lidocaine HCl, Chlorhexidine Gluconate (20% solution), Methyl  
hydroxybenzoate, Propyl hydroxybenzoate, Propylene glycol, Hydroxyethyl  
cellulose, Purified water.

### WARNINGS AND STORAGE CONDITIONS

Keep out of reach of children.  
It must be used only under the supervision of health professional.  
Do not reuse.  
For IM and IV injection.

### Unit In Inner Box Qty in One Box

6ML	25	100
11ML	25	100
12,5GR	20	80