

# ONESEPT PA POWDER HIGH-LEVEL DISINFECTANT FOR MEDICAL INSTRUMENTS AND ENDOSCOPES



EN 13624  
EN 13727  
EN 14476  
EN 14651  
EN 14563  
EN 14562  
EN 14348  
EN 17126  
EN 13704

## USAGE

OneSept PA POWDER is used to provide high level disinfection of all kinds of invasive surgical instruments, dental surgical instruments, surgical laparoscopy instruments, flexible and rigid endoscopes. It is a concentrated powder product.



## PROPERTIES

It is a concentrated powder for professional use  
Anticorrosive effective  
Performs cleaning and disinfection at the same time  
Aluminum, lead, zinc, chrome, copper and nickel/silver planted materials also are not suitable for use  
Solution can be used for 1 business day.

**%99.9**  
Effective

	Test Method	Microorganism	Effect Duration
<b>Bactericidal</b>	EN 13727	Staphylococcus	<b>10 min</b>
		Pseudomonas aeruginosa	
		Enterococcus hirae	
	EN 14651	Staphylococcus	
		Pseudomonas aeruginosa	
<b>Tuberculocidal / Mycobactericidal</b>	EN 14563	Mycobacterium avium	
		Mycobacterium terrae	
	EN 14348	Mycobacterium avium	
<b>Fungicidal / Yeasticidal</b>	EN 14562	Candida albicans	
		Aspergillus niger	
	EN 13624	Candida albicans	
		Aspergillus niger	
<b>Virucidal</b>	EN 14476	Poliiovirus Type 1	
		Human Adenovirus Type 5	
		Murine Norovirus	
<b>Sporicidal</b>	EN 17126	Bacillus subtilis	
	EN 13704	Bacillus subtilis	

## MICROBIOLOGICAL EFFECT

**Bactericidal, Fungicidal, Virucidal and Sporicidal effects.**

## COMPOSITION

Peracetic Acid (solution diluted with %2 water) 2250 ppm,  
Sodium Percarbonate, TAED, Anti-Corrosion, Stabilizer, Auxiliaries.

## WARNINGS AND STORAGE CONDITIONS

Keep out of reach of children.  
Keep away from sunlight and high temperatures.  
Store in tightly closed original packaging/container, dry and cool area.  
Keep away from heat, sparks and fire.  
Under normal conditions, shelf life is 2 years.  
Use opened product within 3 months.  
Store locked up.

KAF G26 - 1	<b>2 KG</b>
KAF G26 - 2	<b>2,5 KG</b>
KAF G26 - 3	<b>5 KG</b>



Peracetic Acid



It is a concentrated powder for professional use



Solution can be used during 1 working day